Harvard PremiumFill[®]

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Harvard PremiumFill[®] is a light cure nano-filled composite for restorations with high demands for esthetics and durability. Beside its low shrinkage and high wear resistance Harvard PremiumFill[®] shows excellent mechanical properties. And above all Harvard PremiumFill[®] has excellent handling properties for a fast and easy application.

Harvard PremiumFill[®] is based on urethanedimethacrylate resin and inorganic filler particles <1.0 µm. The total filler content is 82 % by weight and 74 % by volume.

Harvard PremiumFill[®] is offered in a wide range of Universal, Dentin- and Enamel shades as well as in translucent shades.

Color coded syringes and OptiTips[®] facilitate their classification.

Indications
Suitable for all cavity classes.

Contraindications

The placement of Harvard PremiumFill[®] restorations is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Harvard PremiumFill[®] if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, Harvard PremiumFill[®] may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor allow not polymerized composite to come into contact with materials containing eugenol.

Application

1. Preparation of tooth and color matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate Harvard PremiumFill[®] shade.

Harvard PremiumFill[®] is available in a wide variety of VITA[®]-Shades and in different opacities.

For easy and time saving esthetic restorations the following Universal Shades are available:

A1 U	B1 U
A2 U	B2 U
A3 U	C2 U
A3.5 U	

For highest esthetic demands the shade layering technique is recommended. A wide variety of shades and opacities are here available:

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2. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Harvard CalciumHydroxide).

3. Application of a bonding agent

Apply a common adhesive for enamel and dentine bonding according to the corresponding manufacturer instructions.

The application of the adhesive Harvard Interlock is recommended. Follow the instructions for Harvard InterLock.

4. Application and curing of Harvard PremiumFill[®]

Apply Harvard PremiumFill[®] in the selected shade in increments and light cure. Harvard PremiumFill[®] is easy to model. Place it by using a suitable instrument. Matrix strips may be used.

In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement (e.g. Harvard Ionoresin Fill) or a light cure flowable composite (e.g. Harvard PremiumFlow[®]), followed by the application of Harvard PremiumFill[®].

In case of small cavities Harvard PremiumFill[®] can be applied directly.

For an optimum result apply Harvard PremiumFill[®] in layers of max. 2 mm. Light cure each increment separately.

By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm², cure each increment separately. Hold the light emission window as close as possible to the filling material.

Shades	Incremental depth	Cure time
Universal, Enamel, Translucent	2,0 mm	20 seconds
Dentine	1,5 mm	40 seconds

5. Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs. Polish to high gloss e.g. with polishing discs or polishing brushes.

6. Additional informations

- If Harvard PremiumFill[®] in OptiTips[®] is applied directly into the mouth, for hygienic reasons the OptiTip[®] is for single use only.
- Do not use any resin to adjust the viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard PremiumFill[®] may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA[®] is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!
Store product out of reach of children!

Harvard PremiumFill[®] conforms to the ISO 4049.

The instructions for use has to be kept for the duration of the application.

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Siempre que sea lícita legalmente una exención de responsabilidad, Harvard Dental International GmbH no