CHARISMA® Opal Flow

66076884/11

- (GB)** Instructions for use
- (ES)** Instrucciones de uso
- (BR)** Instruções de uso
- (GR)** Οδηγίες χρήσης
- (CZ)** Návod k použití
- (HU)** Használáti utasítás
- (LV)** Lietošanas instrukcija

- (LT)** Naudojimo instrukcija
- (PL)** Instrukcja obsługi
- (HR)** Upute za uporabu
- (RO)** Instrucțiuni de utilizare
- (TR)** Kullanma talimatı
- (RU)** Инструкция по применению
- (AE)** تعليمات الاستخدام

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Manufacturer:
Kulzer GmbH
Leipziger Straße 2
63450 Hanau (Germany)
Made in Germany

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CHARISMA® Opal Flow

Instructions for use **(GB)**

Product description: Charisma Opal Flow is a flowable, light-curing, radiopaque nano-hybrid composite used for adhesive, tooth-coloured anterior and posterior restorations.

Composition: Charisma Opal Flow is based on multifunctional methacrylate monomers (UDMA, EBADMA) and contains approximately 65% by weight or 41% by volume inorganic fillers such as Ba-Al-F silicate glass, YbF₃, and SiO₂. The filler particle size is between 0.02 µm and 5 µm.

Shades: A1, A2, A3, A4, B2, OM (opaque medium), Baseline

Indications: • Enclosed fissure sealing • Cavity lining - as the first layer for Class I and II cavities • Class V fillings • Minimally invasive Class I and II fillings in areas not subjected to masticatory forces • Minimally invasive Class III fillings • Small repairs of direct and indirect restorations combined with a suitable bonding agent • Splinging of mobile teeth

Contraindications:

Application to dentine close to the pulp (less than 1 mm). Do not use in cases of a known allergy to any of the constituents of Charisma Opal Flow

Side effects:

This product or one of its components may in particular cases cause hypersensitivity reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

Warnings: May cause sensitization with skin contact. May irritate eyes. Avoid contact with the skin and eyes. In case of contact with skin, wash skin with soap and water. In case of contact with eyes, rinse with water and contact a doctor.

The tips supplied for precise, hygienic application are intended for single use and must be discarded after use.

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