HARVARD CEMENT

schnell härtend

0482

HARVARD CEMENT schnell härtend Zinkphosphat-Zement

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HARVARD CEMENT quick setting Zinc phosphate cement

Instruction for use for the dentist Composition:

1 g powder contains:
890 mg zinc oxide,
89 mg magnesium oxide;
1 g liquid contains:
540 mg o-phosphoric acid.

Areas of Application:

Fixation of crowns, bridges and inlays.

Dispensing Instructions:

Fixation consistency:

1.8 g powder : 1.0 g liquid

Dispensing and testing climate for testing in accordance with DIN EN ISO 9917: Mix quantities indicated at 23 ± 1°C and a relative humidity of 50 ± 10%.

Processing Recommendations:

A clean, dry glass plate at ca. 20°C (68 °F) and a non-discolouring spatula are recommended for mixing.

Mixing:

In order to achieve a homogenous consistency, divide the measured quantity of powder into 4 portions (1/2 1/4 1/8 1/8). Beginning with the smallest quantity, stroke the powder into the measured liquid within 90 seconds. The fixation consistency has been reached when the peak formed on lifting the spatula slowly falls back into the paste.

Discard remaining powder or impure liquid.

Mixing Time:

90 seconds

Processing Time from start of mixing

(at 23 °C, 72 °F):

2 minutes 30 seconds

Setting Time from end of mixing:

3:00 – 5:00 minutes

Special Information:

- The mixed cement should never be applied directly onto dentine near to pulp or onto exposed pulp.

- Irritations of the pulp may occur following the cementing.

- The cement must be applied to the dry tooth surface and must set without contact to moisture.

- Close containers immediately after removing material.

- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.

Shades:

1 (white), 2 (bluish white), 3 (yellowish white), 4 (light yellow), 5 (yellow), 8 (pearl grey), 12 (brown), 15 (pink)

Package Sizes:

Powder: 100 g and 35 g

Liquid: 40 ml and 15 ml

HARVARD CEMENT quick setting

- Do not use after the expiration date.

- Storage temperature: 2-28 °C (36-82 °F).

Only for dental use!

Store product out of reach of children!

- Compliant with DIN EN ISO 9917.

The instructions for use has to be kept for the duration of the application.

The user is responsible for correct application of the product according to this instructions for use.

Powder	Order no.	
	Clinic power 100 g	Single powder 35 g
1 - white	7001201	7001501
2 - bluish white	7001202	7001502
3 - yellowish white	7001203	7001503
4 - light yellow	7001204	7001504
5 - yellow	7001205	7001505
8 - pearl grey	7001208	7001508
12 - brown	7001212	7001512
15 - pink	7001215	7001515
	Clinic liquid 40 ml	single liquid 15 ml
Transparent	7001300	7001600

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