**VOCO****Easy Glaze®****GB**

Instructions for use

**CE** 0482**Product description:**

**Easy Glaze** is a nano-filled, light-curing coating with natural fluorescence for surface sealing. The coated surface is lustrous, provides protection against discolouration and yields a natural appearance.

**Easy Glaze** can be polymerised with all dental light-curing devices.

**Indications:**

- Surface sealing of provisional crowns and bridges, glass ionomer restorations and definitive composite restorations
- Sealing and protecting the adhesive interfaces between restoration and tooth structure
- Protecting glass ionomer surfaces against the effects of moisture and dehydration immediately after placement
- Sealing glass ionomer liners/build-up restorations before taking impressions

**Preparation:****Preparation of provisional crowns and bridges**

1. Remove inhibition layer with alcohol
2. Finish and fit the temporary
3. Clean and dry surfaces

**Preparation of glass ionomer materials**

Apply material immediately after finishing the glass ionomer material (see **Application of Easy Glaze**).

**Preparation of a new restoration**

1. Finish the restoration with a superfine diamond bur. Alternatively, polish the surface according to the manufacturer's instructions.
2. Spray away material dust and residues with water. Subsequently dry with oil-free compressed air.

**Preparation of an existing restoration**

1. Slightly roughen the surfaces with a fine diamond bur.
2. Spray away material dust and residues with water. Subsequently dry with oil-free compressed air.

**Preparation of preparation margin**

Uncut enamel should be etched with 35 - 40 % phosphoric acid for 15 s. Subsequently rinse and dry with compressed air.

**Application of Easy Glaze:**

1. Ensure that the surfaces to be treated are dry.
2. Dispense a few drops of **Easy Glaze** into a disposable dish. Replace bottle cap immediately after use.
3. Apply material within 1 min. to guarantee that an oxygen inhibition layer does not form.
4. Apply a thin and air bubble-free layer of the material onto the relevant surfaces with the enclosed brush.
5. Light-cure all treated surfaces with a visible light (> 500 mW/cm<sup>2</sup>, halogen or LED) for 30 s immediately after application. Hold the light source as close as possible to the surfaces being polymerised. Repeat light-curing if the surfaces appear tacky or yellow.
6. Repeat steps 2 - 5 if a second layer appears to be necessary.
7. Where appropriate, check occlusion.

**Information, precautionary measures:**

- Do not use in combination with materials containing eugenol as they can impair the polymerisation or adhesion of **Easy Glaze**.
- **Easy Glaze** is volatile. Only use in well-ventilated areas. Keep bottle out of direct sunlight.
- **Easy Glaze** is flammable. Do not use near open flames or ignition sources.
- **Easy Glaze** contains multifunctional acrylates, methyl methacrylate and BHT. **Easy Glaze** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- In case of contact with the eyes, immediately rinse with water and seek medical advice.
- In case of contact with skin or oral mucosa, immediately remove the material with a cotton pad. Thoroughly rinse affected areas with water afterwards.

**Storage:**

Store at room temperatures of 15 °C to 23 °C. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VOCO GmbH