Fig. 1


Fig. 2



Fig. 3

ACCESSORIES AND SPARE PARTS / ACCESSORI E RICAMBI / ZUBEHÖR UND ERSATZTEILE ACCESSOIRES ET RÉCHANGES / ACCESORIOS Y RECAMBIOS AKCESORIA I CZĘŚCI ZAMIENNE / ΕΞΑΡΤΗΜΑΤΑ ΚΑΙ ΑΝΤΑΛΛΑΚΤΙΚΑ	
REF 520.00	Protection shield / Schermo di protezione / Écran de protection / Pantalla de protección / Escudo de proteção / Osłona zabezpieczająca wzrok / Προστατευτικό πλάνισιον
REF 530.01	Ø 8 mm Optical fibre / Fibra ottica / Fibre optique / Fibra óptica / Fibra ótica / Światłowód / Οπτική ίνα

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MOON/TKD502.00/Rev.07.2025

DESCRIPTION
Curing lamp for professional dental use. Intended use: curing lamp for polymerization of light-curing dental materials (composites, compomers, adhesives, cements) used in conservative dentistry, orthodontics and prosthetics. Intended patient population: any person without restriction concerning age, gender, race or culture. Patient contra-indications: there are no specific patient contraindications or side effects, if the device is used as intended. The light source is provided by a single LED which emits in a narrow band of wavelength, from 430 to 490 nm, with a peak on 460 nm. This band is valid for all materials whose photostarter is made up of camphorquinone (in any case, check whether the material to be cured is sensitive to this wavelength). To ease operation, lamp can be used with various setups: straight, pistol or in any intermediate position (Fig.2). The lamp allows 6 curing cycles, different in terms of emitted light intensity and emission times. Description of the lamp (Fig.1): a) Handpiece b) Swivel end part c) Optical fibre d) Shield to protect eyes e) Control keyboard This medical device meets the requirements of the Regulation (UE) MDR 2017/745 (Class I). Before using the device, read the following warnings carefully: - The device must be used by competent and qualified personnel only; - The device must be used according to these operating instructions; - The device must be used only for the purpose for which it is intended; - The device does not produce electromagnetic interference; - The device is not suitable for use in presence of a flammable anaesthetic mixture with air, oxygen or nitrous oxide; - The device is not liquid-proof (IP20). CONTENTS MOON curing lamp, an optical fibre, a protection shield. The complete SET (REF 502.9x) also includes the silicone hose. CONNECTION Attach the curing lamp to the connector of the supply silicone hoses (REF 330.2x) and screw the corresponding locking nut. Lamp must be connected to a dental unit meeting the requirements of CEI EN 60601-1 standard that provides double insulation from the mains supply.
LED WARNINGS
The LED light is a source of a class 2 according to CEI EN 62471 standard. Such light can damage eyes in case of direct emission without protection. Never use the curing lamp without the special protection shield on and pay attention not to direct the light beam into eyes. Patients who suffer from eye pathologies (such as those who have had cataracts removed or who have trouble with the retina) must be adequately protected during lamp use, for example by wearing appropriate protection glasses. Pay attention to precisely direct the light beam onto the material to be cured. Emitted light can also damage soft tissues (oral mucous membrane, gums, skin). Do not shake the handpiece excessively and pay attention not to drop it, especially the optical fibre. In case of knock or fall, check handpiece integrity before using the curing lamp. Then turn it on and check its correct operation without using it on a patient. In case of cracks, or of any other malfunction, do not use lamp on a patient and contact technical service.
OPTICAL FIBRE
To insert the optical fibre into the handpiece, push it fully down into its housing until a click can be heard. The optical fibre is made in glass and is therefore fragile: in case of knocks or falls it may break or chip, producing sharp bits which may cause injuries and compromising final quantity of emitted light. If the fibre is damaged, a bright dot appears on the point where it is cracked. In all these cases, the optical fibre must be replaced.
DESCRIPTION OF OPERATING CYCLES

CYCLE	POWER EMISSION	Energy	Description / LED (Fig.3)
1	1000 mW/cm² for 20 seconds	20000 mJ/cm²	Standard cycle. LED 1 turns on.
2	1600 mW/cm² for 15 seconds	24000 mJ/cm²	Fast cycle. LED 2 turns on.
3	1800 mW/cm² for 20 seconds	36000 mJ/cm²	Strong cycle. LED 3 turns on.
B	500 mW/cm² for 5 seconds, ramp from 500 to 1000 mW/cm² for 5 seconds and then 1000 mW/cm² for 5 seconds.	11250 mJ/cm²	Bonding. LED 5 + LED 1 turn on. Ramp cycle for a total time of 15 seconds.
R	500 mW/cm² for 5 seconds, ramp from 500 to 2200 mW/cm² for 5 seconds and then 2200 mW/cm² for 5 seconds.	20250 mJ/cm²	Rapid restoration. LED 5 + LED 2 turn on. Ramp cycle for a total time of 15 seconds.
L	500 mW/cm² for 5 seconds, ramp from 500 to 1800 mW/cm² for 5 seconds and then 1800 mW/cm² for 10 seconds.	26250 mJ/cm²	Long restoration. LED 5 + LED 3 turn on. Ramp cycle for a total time of 20 seconds.

Data values are valid as output of the ø8mm optical fibre.

The lamp is equipped with permanent memory so that at power-up the last used cycle is recalled. The lamp is also equipped with a buzzer which emits a beep when light is turned on, every 5 seconds of operation, and when light is turned off.
Furthermore, the lamp is equipped with a thermal protection.

Sequence of operations
Extract lamp from its holder on the dental unit.
Rotate lamp and/or optical fibre to the best setup for use (straight, pistol or intermediate position).
Select the cycle to be used by means of MODE button.
Put optical fibre into the appropriate position: the output of the optical fibre must be kept as close as possible to the material to be cured, though without touching it.
Start the cycle by pressing the START button and wait until light emission turns off automatically.

MODE button (Fig.3)
Press this button to select the cycle to be used. The set cycle is displayed by the turn-on of the corresponding LED (or pair of LEDs). Cycle selection is only possible when lamp is not emitting light.

START button
Press this button to start the set cycle. If this button is pressed again during the cycle, light emission will be immediately stopped.

FAULT SIGNALS
In case of malfunction of the curing lamp, the following signals are shown on the keyboard:
• LED 5 and LED 1 on: Light beam not being emitted by the lamp. Contact Technical Service.
• LED 5 and LED 2 on: Microcontroller for starting the lamp is malfunctioning. Contact Technical Service.
• LED 5 and LED 3 on: Insufficient power supply. Check supply voltage.
• LED 5 and LED 5 blinking: Handpiece thermal shutdown has tripped. These LEDs will continue to blink until the lamp has cooled down sufficiently (approx. 5 minutes) to be used again. If the problem persists, contact Technical Service.

CLEANING AND DISINFECTION
Before sterilizing the optical fibre and protection shield, check that there are no residues of cured materials on them; if necessary, remove residues with alcohol or by using a plastic blade.
Curing lamp can be externally disinfected with appropriate products, as for instance germicidal cleansers based on glutaraldehydes, ethanol or propanol.
To disinfect the handpiece use soft disposable paper, avoiding use of corrosive cleansers and avoiding immersing it in any type of liquid.
To disinfect the handpiece, it is advisable to carry out this operation with the optical fibre inserted. Never use any type of disinfectant on the exposed optical surface of the handpiece when the fibre is removed, as the disinfectant could make this surface become irreversibly cloudy.

STERILIZATION
 Either the optical fibre or the protection shield can be sterilized in a steam autoclave at 134 °C/2 bar (273 °F / 29 psi).
Devices are supplied "not sterile".
The optical fibre can withstand 500 sterilization cycles in autoclave after which it begins to get cloudy and to emit a smaller quantity of light. The protection shield must also be replaced after 500 sterilization cycles.
The handpiece cannot be sterilized in autoclave.

TRANSPORT AND STORAGE CONDITIONS
Temperature: -10 ÷ 60 °C (14 ÷ 140 °F)
Relative humidity: 10 ÷ 85 %
Atmospheric pressure: 50 ÷ 106 KPa

SERVICE
In case of malfunction of the device, and for any overhaul and repair work, please contact your usual supplier or directly to TKD.

REFUSE DISPOSAL
 This device needs to be recycled. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards. The user must return the equipment to its dealer or establish direct contact with an approved body able to process and derive value from this type of equipment (European Directive 2002/96/EEC).

GUARANTEE
The device carries a 24-month guarantee against all defects of construction. Guarantee validity is subject to specific conditions.

INCIDENTS
If any serious incident occurs in relation to the device, report it to a competent national authority, as well as to the manufacturer through your usual distributor.

TECHNICAL SPECIFICATIONS
– Conformity: ISO 10650, CEI EN 60601-1, CEI EN 60601-1-2, MDR 2017/745
– Degree of protection against direct or indirect contacts: class 1, type-B applied part according to CEI EN 60601-1 standard
– LED emission wavelength: 430 ÷ 490 nm (peak at 460 nm)
– LED classification: class 2 source, according to CEI EN 62471 standard
– LED lifespan: 5400000 cycles
– Maximum power emission: 2200 mW/cm² ± 10%
– Maximum curable thickness per single cycle: 3 mm
– Handpiece weight: 122 g
– Supply voltage: 20 ÷ 36 Vdc or 24 ± 10% Vac, 50/60 Hz, Maximum current consumption: 350 mA
– The power supply unit must meet the requirements of CEI EN 60601-1 standard
– Working temperature: 10 ÷ 40 °C (50 ÷ 104 °F)
– Working relative humidity: 30 ÷ 80 %

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