

Specyfikacja Techniczna DMETEC

Specifications of unit:			
Model:	Compact S, Compact S(LED)		
Supply voltage:	AC/DC 24V		
Supply frequency:	50 / 60 Hz		
Power input:	24 VA		
Output characteristics(Ultrasonic output):	Maximum output power with load : 24VA Frequency range available : 24-32kHz.		
Water supply pressure:	1 to 5 bar (15 to 72 psi)		
Dimensions / Weight :			
Circuit module:	Weight:	110g	
	Dimensions:	height:	34mm
		width:	60mm
		depth:	54mm
Handpiece:	Weight:	41g	
	Dimensions:	diameter(max)	17mm
		length	104mm